

DreamWear patient fitting guide

Achieving the right fit with DreamWear is easy. Simply follow the instructions below to get started.

Four cushion sizes (Small, medium, large, medium wide) are available.
Please read instructions for use for further fitting information.

 MultiDoctorShop
Medical products at home

We provide the official document of the manufacturer
www.multidoctorshop.com



1 Hold the sizing gauge horizontally under the nose and select the cushion size based on where the outer edge of the nostrils and tip of the nose contact the gauge.



2 The cushion opening will be directly under the nostrils. The cushion will hug the nose and leak will be minimal.

(At no time should the nose be inside the cushion opening. If any part of the cushion overlaps the nose, the cushion is being worn incorrectly.)




3 With the mask assembled, place the cushion under the nose.



4 Position the frame on the top of the head.



5 Pull the headgear over the back of the head.



6 To adjust the mask, peel the headgear tabs away from the fabric. Adjust the straps length. Press the tabs back against the fabric to reattach.

(Note: Do not overtighten the headgear.)

7



Position the mask until it fits comfortably. When finished, the elbow should rest at the top of the head.

Mask frame sizing

The medium (MED) mask frame will comfortably fit most faces. If the MED frame does not fit a patient's face, a small (SM) or large (LG) mask frame might better suit your patient's needs. **Call for more information, 1-866-992-7279**

Part of the **Dream Family**

DreamWear is part of the Dream Family from Philips Respironics. The Dream Family offers innovative, comprehensive sleep therapy technology like our sleep therapy system, DreamStation, and our patient self-management tool, DreamMapper.



Call for more information, 1-866-992-7279