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Congratulations on your purchase and thank you for choosing our unit.  
 Visit our Internet site [www.flaemnuova.it](http://www.flaemnuova.it) to view the whole range of Flaem products

**Carefully read these instructions and keep them in a safe place for future reference.  
 Use the unit only as described in this instruction manual.**

## IMPORTANT SAFEGUARDS

**This a medical device for home use. It is to be used only upon medical prescription.**  
 Before using the unit for the first time, and periodically throughout its life, check the power supply cord for any damage. If there is any damage, unplug immediately and return the product to an authorized dealer for service.

- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- The unit should not be used in the presence of anaesthetic substances inflammable with air, oxygen or nitrous oxide.
- Never use adaptors for power supplies different from the one indicated on the bottom of the unit. Keep the cord away from heated surfaces.
- Do not handle the plug with wet hands. Never use the unit while having a bath or shower. Never submerge the unit in water; should this accidentally occur, unplug immediately - do not reach for the submerged the unit, always unplug first. After removing a unit from water, do not use again, ensure it is returned to an authorized service dealer.
- Do not wash the unit under running water or do not submerge it in water.
- Unit with case not protected against ingress of liquids.
- Always unplug the unit immediately after use.
- Repairs should be performed by authorized personnel only. Unauthorized repairs will invalidate the warranty and may represent a hazard to the user.

## INSTRUCTIONS FOR USE

**Before each use, the nebulizer and the accessories should be cleaned according to the instructions in the section "CLEANING AND DISINFECTION".**

1. Plug the power cord into a power socket corresponding to the voltage of the unit.
2. Carefully wash your hands before preparing the inhalation therapy.
3. Open the medication nebulizer by turning the upper part (9c) counterclockwise.
4. Pour the drug prescribed by the physician into the lower part (9a). Close the medication nebulizer by turning the upper part (9c) clockwise.
5. Assemble the desired accessories as shown on the drawing (picture 1) (page 1).
6. Sit in a comfortable position while holding the nebulizer, place the mouthpiece in your mouth or use the nasal prong or the mask. Start the unit by turning the switch on (3). Drug delivery to the airways can be maximised by coordinating breathing with the intermittent closure of the drug saving device (8). Close the hole of the drug saving device (8) with your finger and gently inhale as deep as possible. Hold your breath for a moment in order to allow the inhaled aerosol to deposit; then slowly breath out. Whilst holding your breath, remove the finger from the drug saving device (8) to minimise drug wastage.

**The drug saving device (8) should not be used if you require continuous nebulization especially with children or disabled persons.**

7. When the treatment is finished, turn off and unplug the unit. Before reusing, ensure that it has cooled down to room temperature.
8. It may happen that during the use of the nebulizer visible humidity deposits form in the connection tube (6), disconnect the tube from the nebulizer and let it dry using the compressor air flow.

## ACCESSORIES - TECHNICAL FEATURES

### Nebulizer RAPIDFLAEM 4 (suitable for intensive use)

The RAPIDFLAEM 4 is a professional, durable and rapid delivery nebulizer suitable for patients with chronic pathologies especially when costly drugs are used. The constant average size of the aerosol particles delivered by the RAPIDFLAEM 4 nebulizer makes it suitable for treatment of both the upper and the lower respiratory tract. It features a valve system that is automatically controlled by the patient's respiratory rhythm. During inspiration, the inspiratory valve (A) in the nebulizer opens to allow the intake of additional air which increases the amount of nebulized drug available to the patient, while the exhalation valve (B) in the mouthpiece (10) is closed. The patient can keep the mouthpiece in the mouth when breathing out, as the exhalation valve (B) automatically opens when exhaling and the inspiratory valve (A) is closed, which also reduces the drug nebulization and minimises the release of nebulized drug into the environment.

**Essential oils or balsamic substances should not be used with the nebulizer.**

### Technical features

- Made from non toxic materials
- Disinfect by boiling or, when used in clinics or hospitals, autoclave at 121°C max. (1.2 bars) for maximum 25 minutes.
- **Capacity:** max. 7 ml
- \* **Aerosol output:** 0.65 ml/min approx
- \* **Particle:** MMAD 2.45 µm

\*Average values measured with 0.9% saline solution.

Certified with TSI Aerosizer 3220 computerized laser system in compliance with Flaem internal procedure.

## CLEANING AND DISINFECTION

### Accessories

Carefully wash all parts used with the medicament under running water with mild (non abrasive) detergents or in the dishwasher.

All accessories (except for the connecting tube) can be disinfected or sterilised choosing one of these methods:

- a) by boiling in a bain-marie (for max 5 min.); this method cannot be used for the mask
- b) by soaking in a solution of water and vinegar (3/4 water, 1/4 vinegar). Rinse thoroughly with drinking water.
  - The RAPIDFLAEM 4 nebulizer is autoclavable (see 'Technical features')
  - Should the valves in the RAPIDFLAEM 4 nebulizers become blocked, disassemble the individual parts, clean and disinfect them as described above.

### Unit

Wipe the case only with a damp cloth using non abrasive detergents.

**The unit is equipped with an air filter (5) which should be replaced when dirty. Use only Flaem's original accessories.**

## IMPORTANT!!

### Microbial contamination

A personal use of the nebulizer and of the accessories is recommended in order to avoid any risks of infection.

### Interactions

The materials used that come into contact with drugs have been tested with a wide range of drugs. However, due to the variety and continuous evolution of the drugs, the possibility of a chemical interaction cannot be excluded. Once opened, the drug should be used as soon as possible, and prolonged contact between drug and nebulizer should be avoided.

**Disposal:** to safeguard the environment, when the appliance has reached the end of its life take it to an authorized electric appliance disposal center where the raw materials can be recycled.

## UNIT TECHNICAL FEATURES

### Mod. NEBULFLAEM F2000

Voltage:	230V ~ 50Hz 210VA	115V ~ 60Hz 190VA	220V ~ 60Hz 200VA	100V ~ 50/60Hz 200VA
Fuse:	T 2A - 250V	T 3A - 250V	T 2A - 250V	T 5A - 250V
In compliance with:	Dir. 93/42/EEC (EN 60 601-1/A2)			
Approvals:	 0051 			

Operating pressure (with nebulizer):	1.1 bar approx.
Max pressure:	3.5 bar approx.
Compressor air output:	14 LPM approx.
Dimensions:	20x27x10.5 cm
Weight:	2.300 Kg
Sound level (at 1 m):	55 dB (A)
Continuous use	
Operating requirements:	temperature min. 10°C; max 40°C air humidity min. 10%; max 95% <b>RH</b>
Storage requirements:	temperature min. -25°C; max 70°C  air humidity min. 10%; max 95% <b>RH</b>
Operating/storage atmospheric pressure:	min. 690 hPa; max. 1060 hPa <b>P</b>

### SYMBOLS

-  Class II unit
-  Type B unit
-  Warning: follow the "Instructions for use"
-  Alternating current
-  Switch OFF
-  Switch ON
-  Never use the unit while having a bath or shower
-  Latex Free

### Electromagnetic compatibility

This aerosol therapy unit is designed to meet the current requirements for electromagnetic compatibility. However if you suspect that the unit is interfering with the operation of your TV, radio or other electrical appliance, try relocating the unit until the interference ceases. Alternatively try plugging the unit into a different power socket.

**Warranty:** The unit is under warranty from the dealer that sold you the unit.

*The Manufacturer, the Retailer and the Importer are responsible for safety, reliability and performance of the unit only if:* • the unit is used in accordance with the instructions for use; • the electrical system of the facility where the unit is used adheres to IEC standards; • assembly, additions, recalibrations, modifications or repairs are carried out by authorized personnel.

The manufacturer reserves the right to introduce technical changes.

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© Copyright by FLAEM NUOVA S.p.A., 2004