



NOBA Verbandmittel

Danz GmbH u. Co KG

NOBATON®

Product Description and Purpose

Finely woven, tightly knitted tubular bandage, which is made of unbleached cotton. A high crosswise elasticity ensures a good adaptability. The unbleached cotton shows water-repellent features and convinces by its excellent skin tolerance. NOBATON® serves as underlay under rigid or half-rigid cast bandages. Its tight weave makes it especially suited for usage with the half-rigid cast bandage NOBASOFT®. The cast bandage NOBASOFT® can be applied directly on the tubular bandage whose tight structure protects the skin against the material soaked with polyurethane resin. It also allows the application of a tightly fitting bandage.

Composition

Cotton

Normative and Legal Requirements

NOBATON® is a medical product in accordance with the MDD, class I, rule 1, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and EN 980 labels on all its packaging.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

