



NOBA[®] - EXTENSIONSVERBAND

with adhesive

Product Description and Purpose

The ready-to-use skin traction kit is provided with two strongly adhesive plaster strips. It serves as an extension in the case of fractures of the thigh or lower leg bones. The skin traction kit is adhered to the lower leg and additionally fixed with the enclosed bandage. Weights can be attached to the affixed strings via a suitable roller system. The bandage is adequately padded in the ankle region.

Composition

Polyester, cotton, polyurethane, polypropylene, polyamide, synthetic zinc oxide adhesive

Normative and Legal Requirements

NOBA[®] - EXTENSIONSVERBAND is a medical product in accordance with the MDD, class I, rule 1, and complies with the requirements of the directive 93/42/EEC.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: bandages in foil made of cellophane and folding box made of cellulose

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

