

## Mirage Vista™ Nasal Mask



### The unobtrusive choice

ResMed's Mirage Vista™ Nasal Mask combines the proven comfort and performance of ResMed masks with an unobtrusive, streamlined design.

The unique design increases field of vision, and allows patients to wear glasses if required.

### User-friendly

A number of innovative features make the Mirage Vista easy to fit, put on and take off, disassemble, reassemble, and clean.

### Performance and comfort

The Mirage Vista uses the proven and popular ResMed Mirage™ cushion technology. This, together with unique headgear, ensures stability and seal with long-term comfort.

Mirage Vista Nasal Mask



**Stabilizing panels**

hold headgear in correct position

**Headgear**

has easy-to-use headgear clips and a non-slip crossover buckle

**Elbow clip**

provides convenience & ease of use

**Angle selector**

provides fine tuning for optimal fit

**Easy-lock parts**

allow easier handling

### Mirage cushion technology

- Trusted performance of ResMed Mirage cushion technology.
- Soft, dual-wall cushion reduces pressure on the nasal bridge.
- Optimizes seal even when patient moves during sleep.



Mirage Cushion Technology—soft but with optimum seal

### Headgear

- No further readjustments necessary after first fitting.
- Stabilizing panels hold headgear in place on the head for added comfort and stability.
- Comes fully assembled for convenience.
- Breath-O-Prene® fabric allows air to flow through and moisture to escape.



Headgear—no readjustment needed

### Elbow clip

- Simply snaps into position for convenience and ease of use.
- Allows the patient to breathe normally while you fit the patient's mask, prior to connecting the air tube.
- Can be easily removed while the mask is still on the face.
- Has 360 degrees of rotation—combined with hose swivel gives increased mobility during sleep.



Elbow clip—easy and convenient

### Easy-lock parts

- More preassembled parts and fewer parts make mask handling easier—even for patients with dexterity problems.



Easy-lock parts—easier handling

## Mirage Vista Nasal Mask System

### Product Codes

#### North & Latin America

Standard **60000**

Deep **60001**

#### Europe & Asia Pacific

Standard **60003**

Deep **60004**



The Mirage Vista Nasal Mask allows patients to follow their normal bedtime routines prior to sleep

### Suitable for a range of patients

Recommended first choice for first-time CPAP users, and particularly suitable for a range of “difficult-to-fit” patients:

- anxious/claustrophobic patients
- asymmetrical facial structures
- patients with dexterity problems
- vision-impaired patients.

### Technology highlight

#### Angle selector

- Offers five optional mask positions for optimum mask comfort.
- Ensures Mirage Vista Nasal Mask conforms to patient's face.
- Left and right sides can be adjusted independently for customization to asymmetrical facial structures.
- Helps provide further stability and seal without need to tighten headgear.

