



## RUDABLOCK<sup>®</sup>-transparent

### Product Description and Purpose

RUDABLOCK<sup>®</sup>-transparent is an individually packed sterile pressure plaster with special wound pad. The wound pad is overlapped on all sides by the carrier material and is thus protected. The pressure plaster is specially designed to staunch blood. If the wound contact layer gets in touch with blood it swells up and exerts pressure on the wound. This pressure expedites the staunching of blood. The special plaster is suited for the staunching of blood in dialysis as well as for smaller wounds (for example cuts, stab wounds).

### Composition

Viscose, polyurethane, polyacrylate adhesive

### Normative and Legal Requirements

RUDABLOCK<sup>®</sup>-transparent is a medical product in accordance with the MDD, class Is, rule 4, and complies with the requirements of the directive 93/42/EEC.

Sterilization of the product complies with DIN EN 11137.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

### Packaging

|                      |                               |
|----------------------|-------------------------------|
| Primary packaging:   | paper-paper packaging         |
| Secondary packaging: | folding box made of cellulose |
| Tertiary packaging:  | carton made of cellulose      |

### Storage

To be stored in a dry and dust-free environment, protected against direct sunlight

The product is bearing the CE marking, in compliance with the directive's provisions

