

# D100

CPAP MASKS



We provide the official  
document of the manufacturer  
[www.multidoctorshop.com](http://www.multidoctorshop.com)

## Nasal and Full Face Mask Solutions



# D100

## CPAP MASKS

The D100 series CPAP masks offer **comfort, durability and simplicity for effective CPAP therapy whilst remaining cost-effective for both patients and providers alike.**

### Features and Benefits

- Available as a nasal mask (D100N) or a full face mask (D100F).
- Choose from three sizes of silicone cushion; S, M and L.
- Three sizes of forehead pad included with every mask for maximum comfort.
- Fully adjustable headgear.
- 360° swivel for optimum tubing positioning.
- Easy to clean and maintain.
- Latex free.

Available in both a nasal and full face design, the D100 masks feature a moulded silicone cushion and are lightweight, durable and comfortable to wear.

The D100N nasal mask and D100F full face mask are available with three sizes of cushion to fit most profiles. The D100 series headgear makes the mask simple to fit and can be easily adjusted to achieve an effective seal, leading to improved patient therapy adherence.

The tubing attachment has a 360° swivel which allows the user to adjust the tubing according to their preferred sleeping position.

The D100 series CPAP masks from DeVilbiss Healthcare provide users with a full range of features and exceptional value.



### SPECIFICATIONS

#### Weight D100N Nasal Mask

with Headgear	88g
Headgear	31g

#### Weight D100F Full Face Mask

with Headgear	122g
Headgear	31g

#### Material

Cushion	Silicone
Forehead Pad	Silicone

#### Storage

Temperature	-18° C to +50° C
Relative Humidity	< 85 %

#### Resistance D100N Nasal Mask

50 LPM	0.47 cmH <sub>2</sub> O
100 LPM	2.24 cmH <sub>2</sub> O

#### Resistance D100F Full Face Mask

50 LPM	0.24 cmH <sub>2</sub> O
100 LPM	0.62 cmH <sub>2</sub> O



D100 Nasal Mask



D100 Full Face Mask



Three sizes of cushion available



D100 Headgear

Distributor:

**DeVilbiss**  
HEALTHCARE

DeVilbiss Healthcare Ltd.  
Unit 3, Bloomfield Park, Bloomfield Road  
Tipton, West Midlands DY4 9AP  
United Kingdom  
Phone + 44 (0) 121 521 3140  
Fax + 44 (0) 121 521 3141  
www.devilbisshc.com  
enquiries@devilbisshc.com