



**PHILIPS**  
RESPIRONICS

*DreamStation*

CPAP

## DreamStation helps support long-term patient use while creating **new efficiencies for your business.**

**Connectivity** options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth® on every device, and offering optional Wi-Fi and cellular modems.\*

**Performance Check** simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

**Daily Progress Feedback** displays a simple trend of patients' nightly use to help track progress on the path to adherence.

### **DreamStation features**

- Mask type resistance control
- Mask fit check
- Oximetry capable
- Encore patient management system with optional Adherence Profiler
- DreamMapper patient self-management system compatibility

### **Humidification**

- Fixed, adaptive, heated tube
- Pre-heat humidifier option
- Modular design with one-piece, dishwasher-safe water chamber

### **Protection programs**

- Flex Promise program
- Bi-level Rescue program
- "One-Touch" warranty programs

### **CPAP features**

- C-Flex
- SmartRamp
- Basic AHI/Leak reporting
- Standard integrated Bluetooth

 **MultiDoctorShop**  
Medical products at home

We provide the official  
document of the manufacturer

[www.multidoctorshop.com](http://www.multidoctorshop.com)

\*Internal assessment of 2015 competitive CPAP data comparing to ResMed Aircurve10/Aircurve10 platform and Fisher & Paykel Icon series platform. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respiroics is under license. Other trademarks and trade names are those of their respective owners.

# DreamStation CPAP

| Comfort features |          |              | Encore data capture and reporting |                |
|------------------|----------|--------------|-----------------------------------|----------------|
| Modes            | Features | Flex Comfort | Advanced Event Detection          | Flow Waveforms |
| Fixed CPAP       | N/A      | C-Flex       | N/A                               | N/A            |

## Specifications

### DreamStation CPAP

|                                 |  |
|---------------------------------|--|
| Pressure range                  | 4 to 20 cm H2O   |
| Flex pressure relief            | 0 to 3   |
| Ramp time                       | 0 to 45 min (5-minute increments)  |
| Starting ramp pressure          | 4 cm H2O to CPAP Pressure (CPAP)   |
| Humidification                  | Heated humidification: Fixed, Adaptive, Heated Tube modes; Heated Tube (optional)      |
| Data storage capacity (minimum) | SD card: > 1 year On-board: <6 months  |
| Altitude compensation           | Automatic  |
| Filters                         | Reusable pollen; disposable ultra-fine   |
| Device controls                 | LCD, control dial/push button  |
| Warranty                        | 2 years (US)   |
| Dimensions                      | 15.7x19.3x8.4cm (w/o humidifier)<br>29.7x19.3x8.4cm (w/ humidifier)                    |
| Weight                          | 1.33kg/2.94lbs (w/o humidifier)<br>1.98kg/4.37lbs (w/humidifier)<br>Incl. power supply |
| Electrical requirements         | 100 – 240 VAC, 50/60 Hz, 2.0-1.0 A   |

## Ordering information

| DreamStation CPAP                                 | Part number |
|---|-------------|
| DreamStation CPAP                                 | DSX200S11   |
| DreamStation CPAP with humidifier                 | DSX200H11   |
| DreamStation CPAP with humidifier and heated tube | DSX200T11   |
| Power supply                                      | Part number |
| DreamStation 80-watt power supply                 | 1118499     |
| Power cord  | 1005894     |

## Ordering information

| Humidification                             | Part number |
|--|-------------|
| DreamStation humidifier                    | DSXH        |
| DreamStation water chamber                 | 1122520     |
| RP-DS Humidifier Dry Box Assembly          | 1120668     |
| RP-DS Humidifier Flip Lid Seal             | 1120617     |
| RP-DS Humidifier Dry Box Inlet Seal        | 1120613     |
| Tubes                                      | Part number |
| 15MM Std Tube – DreamStation – RP          | PR15        |
| 15MM Heated Tube                           | HT15        |
| Filters                                    | Part number |
| Pollen filter, reusable (1 per pack)       | 1122446     |
| Ultra-fine filter, disposable (1 per pack) | 1122447     |
| Ultra-fine filter, disposable (2 per pack) | 1122518     |
| Ultra-fine filter, disposable (6 per pack) | 1122519     |
| Accessories                                | Part number |
| Link module                                | 1120293     |
| Nonin SpO2, assembly                       | 1121694     |
| Reuseable finger sensor (clip)             | 936         |
| Oximetry finger sensor, adult, flex        | 953         |
| Flexiwrap, sens tape, adult, 25/pk         | 954A        |
| Reuseable finger sensor (soft), medium     | 1062050     |
| Reuseable finger sensor (soft), small      | 1062051     |
| Modems                                     | Part number |
| Cellular modem, oximetry capable           | 100600C     |
| Cellular modem, non-oximetry               | 100605C     |
| Wi-Fi accessory, oximetry capable          | 100700W     |
| Wi-Fi accessory, non-oximetry              | 100705W     |

© 2015 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

www.philips.com



Caution: US federal law restricts these devices to sale by or on the order of a physician.

edoc EB 11/12/15 MCI 4107064  
1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 • 724-387-4000