# Mask solutions to Simplify SleepLife















# Simplifying SleepLife

We sleep for a third of our lives. During sleep we grow, we repair, memories are made; our bodies and minds are rejuvenated.

Sleep feeds our quality of life. Just as our work-life, home-life, love-life and social-life are important to our well-being, our SleepLife is crucial to our health and happiness.

Fisher & Paykel Healthcare is a leading innovator that excels in the treatment of Obstructive Sleep Apnea (OSA).

Our high-performance family of masks simplifies the daily challenges of therapy. They are designed and renowned for: **comfort, seal** and **easy use**.





Fisher & Paykel Healthcare, working with you to navigate patients towards adherence.

# **OUR MASK FAMILY**

The Fisher & Paykel Healthcare range of mask solutions has been designed with innovative features to Simplify SleepLife and support the patient's journey to therapy adherence.



Facial comfort is achieved through unique cushioning technologies which contour the face. Together with innovative headgear designs, greater comfort is today being achieved by lighter masks, softer seals and small-frame profiles.



An optimal mask seal is offered and achieved in each delivery category. All are designed to minimize leakage and maximize therapy adherence.



With a focus on Simplifying SleepLife, our family of masks is designed to be easy in every respect: set-up, cleaning, adjusting and sleeping.

Just fit and forget.

# Fit-and-forget technology

Using our anthropometric database with thousands of facial contours, our designers have created a combination of technologies that adjust to the user. Personalized comfort is achieved through the mixture of these intuitive technologies which require very little manual adjustment. Fit-and-forget means just that – quick and hassle-free fitting, even when it comes to cleaning and reassembly.

# **Out-of-box solutions**

All masks come assembled and ready to use. Our range of mask sizes and designs makes fitting easy. There is no requirement for difficult manual adjustments once the masks have been fitted.





# A MASK IN EVERY DELIVERY CATEGORY


Fisher & Paykel Healthcare's masks are designed to work in harmony with the way patients naturally breathe while they sleep. Some patients are nose breathers, others breathe through their mouth and some are a combination. Patient adherence is better supported by choosing a mask based on how a patient breathes. Only Fisher & Paykel Healthcare offers a mask in all four delivery categories. The diagram below guides you to the right F&P mask for your patient healthcare.

# **DELIVERY CATEGORIES MASK CHOICE DIAGRAM**

# **BREATHING MODE**

# STEP 1 How do you predominantly breathe while you sleep? NOSE MOUTH STEP 2 Are you sensitive to enclosed spaces? NO YES NO

# **DELIVERY CATEGORY**



You naturally breathe through your nose or mouth?

You have a nasal obstruction?

You scuba-dive or snorkel?

- You wear glasses and enjoy reading or watching TV before falling asleep?
- You would prefer the freedom of little or no headgear?
- You are comfortable with the feel of a hat or cap?

- A full face mask will allow mixed delivery as your breathing patterns change through the night.
- An oral or full face mask allows oral breathing.
- If so, a mask that provides a free field of vision may suit you better.
- • If so, a nasal pillows or oral mask may suit you best.
- Nasal or full face masks with headgear are not dissimilar in feel to a hat. Each is designed for lightweight, personalized comfort.
- An oral mask slips into your mouth in a similar manner to a diver's mouthpiece.

HINTS ON MASK CHOICE

# **INFORMATION ON MASK MODELS**

ORAL MASKS Oracle™

NASAL PILLOWS F&P Pilairo™ Opus™ 360 NASAL \*
F&P Eson™
F&P Zest™Q
FlexiFit™

FULL FACE MASK
F&P Simplus™
Forma™



## THE ORACLE

is unique – it delivers pressure exclusively through the mouth with minimal headgear.



# OPUS 360

are small and light masks, delivering pressure directly to the nasal passages.



## ESON, ZEST Q AND FI FXIFIT


nasal masks offer over-the-nose solutions that seal well and are easy to fit.



## SIMPLUS, FORMA AND FLEXIFIT 431

are for those who prefer a mix of nasal and oral breathing.

# Performs - in Tune with you



**F&P Eson** maximizes no-fuss comfort and simplicity by minimizing manual adjustments, parts and noise.

The F&P Eson's three simple components, the RollFit™ Seal, ErgoFit Headgear and Easy Frame, work in harmony to deliver the comfort, seal and easy use that Fisher & Paykel Healthcare masks are known for. Like the F&P Eson, the F&P Zest Q and F&P Lady Zest Q uses an Advanced Air Diffuser. This refinement ensures minimal noise and efficient exhaust dispersal.

# Nasal Pillows Masks

# Sleep Freely on Pillows of Air



The **F&P Pilairo** Nasal Pillows Mask takes its design inspiration from the world's lightest little bird... 'the hummingbird'.

The AirPillow™ Seal self-inflates to gently envelop the nose in a "hover-like way" while providing freedom of movement in conjunction with the new Stretchwise™ Headgear. At only 1.83 ounces (52 grams), it's the lightest little performer in the F&P Nasal Pillows class with the least number of parts.

The F&P Opus is our original mask in this category. The F&P Opus is an alternative option with an innovative 360-degree pivoting ball-and-socket elbow which provides for versatility and freedom of movement.

# Simply Fits + Performs



When engineering the **F&P Simplus**, our designers set out to create a mask that revolutionized Full Face comfort, seal and easy use.

The F&P Simplus incorporates three key components, the RollFit<sup>™</sup> Seal, ErgoForm<sup>™</sup> Headgear and Easy Frame, all designed to work in harmony.

As the name suggests, the RollFit Seal "rolls" back and forth on the bridge of the nose to adjust automatically to reduce pressure on the bridge of the nose. The one frame fits all three sizes and has an Easy-Clip frame attachment. The ErgoForm Headgear is ergonomically designed to self-locate high on the rear of the head, allowing for unrestricted head movement for optimal performance.

The F&P FlexiFit<sup>™</sup> 431 Full Face Mask is our original mask with the Under-Chin design.

Oral Masks

# A Unique Solution for Mouth-Breathers



The unique **F&P Oracle** Oral Mask delivers pressure through the mouth rather than the nose.

Clinically proven to function as effectively as a nasal mask, the F&P Oracle is a perfect solution for mouth-breathers and patients with chronic nasal obstruction.

The F&P Oracle is easy to use and easy to seal. With nothing covering the field of vision, the F&P Oracle provides freedom of movement. It is an Out-of-Box Solution, fully assembled and ready for use.

# For more information please contact your local Fisher & Paykel Healthcare representative

# Manufacturer

Fisher & Paykel Healthcare Ltd 15 Maurice Paykel Place East Tamaki, Auckland 2013

PO Box 14 348 Panmure Auckland 1741 New Zealand

Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com

# Australia

Fisher & Paykel Healthcare Pty Limited 36–40 New Street, PO Box 167 Ringwood, Melbourne Victoria 3134, Australia

Tel: +61 3 9879 5022 Fax: +61 3 9879 5232

#### Austria

Tel: 0800 29 31 23 Fax: 0800 29 31 22

#### Benelux

Tel: +31 40 216 3555 Fax: +31 40 216 3554

#### China

Tel: +86 20 3205 3486 Fax: +86 20 3205 2132

# France

Tel: +33 1 6446 5201 Fax: +33 1 6446 5221

## Germany

Tel: +49 7181 98599 0 Fax: +49 7181 98599 66

## India

Tel: +91 80 4284 4000 Fax: +91 80 4123 6044

# Irish Republic

Tel: 1800 409 011

## Italy

Tel: +39 06 7839 2939 Fax: +39 06 7814 7709

## lapan

Tel: +81 3 3661 7205 Fax: +81 3 3661 7206

#### Northern Ireland

Tel: 0800 132 189

#### Spair

Tel: +34 902 013 346 Fax: +34 902 013 379

#### Sweden

Tel: +46 8 564 76 680 Fax: +46 8 36 63 10

#### Switzerland

Tel: 0800 83 47 63 Fax: 0800 83 47 54

#### Taiwan

Tel: +886 2 8751 1739 Fax: +886 2 8751 5625

# Turkey

Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, P.O. Box 06371 Ostim, Alinteri Bulvari 1161/1 Sokak No. 12-14, Ankara, Turkey

Tel: +90 312 354 34 12 Fax: +90 312 354 31 01

# UK EC REP

Fisher & Paykel Healthcare Ltd Unit 16, Cordwallis Park Clivemont Road, Maidenhead Berkshire SL6 7BU, UK

Tel: +44 1628 626 136 Fax: +44 1628 626 146

# USA/Canada

Tel: +1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

