


Eclipse™ 3 with autoSAT™

The Eclipse 3 is a single solution with both continuous and pulse dose options. It's powerful enough to be your stationary system and mobile enough to ambulate.

The autoSAT® feature of the Eclipse 3 adjusts to your respiratory rate as it increases and decreases. autoSAT provides you unparalleled performance without limiting available oxygen, enabling the unit to automatically adjust to increased oxygen demands that occur as a part of everyday life.



The Eclipse 3 can be your 24/7 all-in-one solution, whether for day or nocturnal needs. Providing up to 3 LPM continuous flow and the largest pulse dose size on the market at 192 mL, the Eclipse 3 can adapt as your oxygen needs change.

- Sleek design is small, lightweight and easy to maneuver
- Telescoping-handled cart facilitates movement in or outside the home
- AC, DC and rechargeable battery enables longer excursion times and extended overnight travel anywhere in the world
- FAA-approved for on-board, in-flight use

UltraSense™ technology

FreeStyle, FreeStyle 5, and Focus are equipped with UltraSense™ technology, which provides a more sensitive pulse dose delivery than competitor devices, even those that include a "sleep mode" setting. The advanced intelligence built into the microprocessor of the Focus ensures fast, reliable pulse dose delivery in the first phase of inhalation, when the highest levels of gas exchange are occurring in the lungs. This gives providers and patients alike the peace-of-mind that with UltraSense, patients are reliably receiving the oxygen they need, when they need it most.

FreeStyle™ 5 with UltraSense™

The lightweight and compact portable, FreeStyle 5 offers five pulse setting options. The pulse-flow delivery it offers enables higher-flow patients to be as active as possible. This model provides you with all the advantages of the standard FreeStyle unit, including proven reliability and durability in a slightly larger platform for greater capacity.



FreeStyle 5 can be conveniently operated and recharged through the AC or separate DC power supply. This allows you to depend

on the flexibility of what FreeStyle 5 offers, whenever and wherever you use the oxygen unit.

The AirBelt™ accessory (which is also compatible with FreeStyle and Focus) is designed to be conveniently worn around the waist and it ensures that extended battery time between recharges is possible. This option maximizes duration time up to 7 hours between battery recharges for those occasions when AC or DC power is not readily available.

- Higher capacity POC at 6.2 lbs (2.8 kg)
- Carrying bag can be converted to a backpack using the optional harness accessory
- Operates on rechargeable battery, AC or DC power
- FAA-approved for on-board, in-flight use

FreeStyle™ with UltraSense™

FreeStyle is the lightest POC up to a 3-setting and offers a uniquely efficient, ultra-compact, and flexible pulse-flow delivery system.

You have the ability to move effortlessly between power sources and no longer worry about running out of oxygen. The universal AC/DC power supply simplifies



the process of operating the unit and recharging the battery. With the internal battery, it is possible to reach up to ten hours of ambulation between recharges!

For patients who want an alternative carrying method for FreeStyle 5 or FreeStyle, an optional backpack harness is available. A carry-all bag for unit accessories is included.

- Operates on AC, DC, and rechargeable internal battery power
- Optional AirBelt extends duration up to 10 hours between recharges
- Carrying bag can be converted to a backpack using the optional harness accessory
- FAA-approved for on-board, in-flight use


Focus™ with UltraSense™

At less than 2 lbs (0.8 kg), O₂ patients like you can be more mobile without having to carry burdensome extra system weight. The pulse flow setting of 2 gives you the ability to ambulate like never before!

Simply plug the unit into an AC or DC power outlet, or operate it on battery for unlimited oxygen. Vacations, trips, and other special events are now possible — whether just around the corner or around the world.

Focus is the most miniaturized and convenient medical device for reliable oxygen-making capability. This POC is also packaged with a shoulder strap so that the unit may be easily worn around the waist or carried, depending on preference. With the optional AirBelt, it can extend battery time an additional four hours between recharges.

- Ultra wearable and weighs only 1.75 lb (0.8 kg)
- Operates on external micro-batteries and can be used with AirBelt for extended battery duration
- FAA-approved for on-board, in-flight use



Accessories equipped with
Eclipse 3, FreeStyle 5, FreeStyle, and Focus:

Eclipse 3

- AC and DC Power Supply
- Accessory Bag
- Power Cartridge
- Universal Cart

Optional Accessories:

- Desktop Charger
- Eclipse 2 Plus Firmware Upgrade
- Eclipse Protective Cover
- Eclipse Data Acquisition Tool (EDAT)
- Eclipse Wheelchair Pack
- Fingertip Pulse Oximeter
- Humidifier Adapter Kit
- Travel Accessory Kit
- Travel Case
- Universal Power Plug Set

FreeStyle 5

- AC Power Supply with Power Cord
- Carry-All-Bag for Accessories
- Carrying Case and Shoulder Strap
- DC Power Supply

Optional Accessories:

- AirBelt and AC Power Supply with Power Cord
- Backpack Harness

FreeStyle

- Carry-All-Bag for Accessories
- Carrying Case and Shoulder Strap
- DC Power Cable
- Universal AC/DC Power Supply with Power Cord

Optional Accessories:

- AirBelt and AC Power Supply with Power Cord
- Backpack Harness

Focus

- AC and DC Power Cords
- Carry-All-Bag for Accessories
- Carrying Bag with Battery Pockets
- Carrying Case and Shoulder Strap
- Coil Cord which Connects Battery to Power Supply Charger to Battery
- Coil Cord with Switch that Connects Battery to Focus
- Lithium ion rechargeable micro-battery (2) and Battery Case (1)
- Universal AC/DC Power Supply with Power Cord

Optional Accessories:

- AirBelt and Battery Belt Interface Cable

Eclipse 3, FreeStyle 5, FreeStyle, and Focus may be offered for purchase, as well as through short- or long-term rentals by your local oxygen provider. These products are covered under a variety of payment or insurance options with full or partial reimbursement. Find out if these oxygen options are the right solution for you or a loved one by seeking a demonstration from your local provider.

Regain independence and feel the difference a lighter, O₂ system can make for ambulation.



CAIRE, Inc.	Chart BioMedical, Ltd.
2200 Airport Industrial Way,	Unit 6, Ashville Way
Suite 500 Ball Ground, GA 30107	Wokingham,
800.482.2473 Telephone (CAIRE)	Berkshire, RG41 2PL
800.874.0202 Telephone (AirSep)	Ph +44(0) 1189 367060
888-932-2473 Fax (CAIRE)	Fax +44(0) 1189 799245
716-691-4141 Fax (AirSep)	

Copyright © 2013 Chart Industries
Chart Industries reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or, processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.



The Smallest and Lightest Oxygen Systems — from Palm Size Portables to Continuous Flow Transportable Oxygen Concentrators



Eclipse™ 3 with autoSAT™
FreeStyle™ 5 with UltraSense™
FreeStyle™ with UltraSense™
Focus™ with UltraSense™



THE ART OF
OXYGEN

Eclipse™ 3 Transportable Oxygen Concentrators
FreeStyle™ 5 • FreeStyle™ • Focus™ Portable Oxygen Concentrators