

## Designed **with care**



As the leader in mask technology, ResMed has partnered with pediatric health experts to bring you Pixi—a pediatric mask designed with children and their families in mind.

### A mask designed **specifically for children**

Uniquely constructed using the anthropometric data of children aged two and over, Pixi Pediatric Mask is not merely a scaled down adult mask.

Taking into account the bone structure, facial characteristics and skin sensitivity of children, the mask is designed to minimise pressure on soft faces.

Pixi Pediatric Mask is made with both the child and their carer in mind.

Every feature is carefully considered to address and improve therapy **comfort, acceptance and overall experience for everyone involved.**

Confident that their child is being treated effectively, carers can be free to focus on emotionally supporting their child on their journey.



**Pixi™**  
PEDIATRIC MASK



## Technology Highlight

### NEW Headgear designed especially for children\*

**Three points of adjustment**—the mask can be adjusted at the top, side and bottom straps, providing a greater range of adjustment to cater for children of various ages.

**Out-of-sight**—the headgear is designed to sit away from the child's eyes and ears to reduce irritation and obstruction.

**Clear line of sight**—not only can the child see clearly while wearing the mask, but carers are also reassured and comforted by being able to see their child's eyes and more of their face.

*\*Designed using the specific head and facial anthropometric data of children.*

### Optional pressure port with tubing

not attached to the mask, packed separately

### Soft, comfortable cushion

single-wall, thin silicone cushion minimises pressure on children's soft faces, without compromising seal

### Two tubing positions

air tubing can be moved to either side of the cushion depending on child's sleeping position

### Spring-flex lightweight tube

stretches and bends easily, so it won't get caught or drag as the child moves around during sleep



### Comfortable Breathoprene™ headgear materials

lightly coloured for a more reassuring, unobtrusive appearance on the child's face

### Side-facing back strap

ergonomically designed from the fitting adult's point of view

### Flexible stabilisers

allow the mask to conform to each child's unique facial structure and shape for the most customised fit

### Quick-release latch

makes it easy to take the mask off, while retaining desired headgear settings

### Emergency latch

releases the latch even faster in an emergency, by pulling firmly on the lower headgear strap

## Complete System

|      | Europe 1 and 2 | Europe 3     | USA, Asia Pacific, Australia and New Zealand |
|------|----------------|--------------|--|
| Pixi | <b>61031</b>   | <b>61032</b> | <b>61030</b>                                 |

**Europe 1:** English (UK), French, German and Italian (Switzerland)  
**Europe 2:** English, Danish, Finnish, Norwegian and Swedish  
**Europe 3:** English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish



**Pixi™**  
PEDIATRIC MASK