



# Therapy you can trust

## Philips Respironics Porta-neb nebulizer compressor system

The Philips Respironics Porta-neb is a reliable and easy to use nebulizer compressor system, suitable for nebulizing most aerosol medications prescribed for common respiratory diseases.

For many years, Porta-neb has provided reliability, fast treatments and high quality aerosol therapy for patients.

Coupled with the SideStream Reusable handset, Porta-neb provides a highly efficient aerosol therapy system. This combination results in fast treatment times and a high quality aerosol that will deliver 80% of particles below 5 microns\*, allowing for more drug to reach the target areas of the lungs.

Suitable for professional or home use, Porta-neb has been designed to easily manage prolonged use in a clinical environment as well as frequent use in the home.

### Key features

- Matched with SideStream for enhanced performance
- 7 l/min. high flow compressor
- Continuously rated for constant use
- Contoured design for easy cleaning
- Large integral storage
- 3 year warranty

**PHILIPS**

**RESPIRONICS**

sense and simplicity



Product specification	
Weight	2.9 kg
Dimensions (cm)	22 x 38 x 12
Electrical requirements	230V 50hz, 210 VA
Maximum nebulizer capacity	10 ml
Noise level	55 dBA
Treatment time (2.5 ml charge)	4.2 mins
Dynamic flow	7 lpm
Max pressure	3.5 bar
Operating pressure	1.05 bar
Neb rate	0.41 ml/min
MMD	3.06 µm
Warranty	3 years

The Porta-neb has instructions for use in the following languages: English, French, Italian, German, Portuguese, Spanish, Norwegian, Dutch, Greek, Turkish, Czech, Polish, Swedish, Danish, Finnish and Arabic.

## Reputation for delivering trust

Porta-neb delivers excellent performance and a dependability that professionals trust to help improve the quality of the lives of those affected by respiratory diseases.

### Porta-neb part numbers

Product description & literature languages	Part number
Porta-neb SideStream compressor system (Euro power cord, EN and AR)	1800
Porta-neb SideStream compressor system (UK power cord, EN and AR)	1802 and 1837
Porta-neb SideStream compressor system (Euro power cord, EN, FR, DE, NL, PL, CZ and IT)	1805
Porta-neb SideStream compressor system (Euro power cord, NO and SE)	1806N
Porta-neb SideStream compressor system (Euro power cord, FI and DK)	1807
Porta-neb SideStream compressor system (Euro power cord, EN, ES, PT, TR, PL, CZ and EL)	1811
Porta-neb SideStream compressor system (Australian power cord, English language)	1861

\* Data held on file at Philips Respironics  
 Respironics and Porta-neb are trademarks of Respironics, Inc., all rights reserved.

### Porta-neb spares and accessories

Product description	Part number
Reusable SideStream handset kit containing SideStream nebulizer, adult mask, pediatric mask, reusable mouthpiece and supply tubing	2010A
Supply tubing	2150A
Supply tubing (pack/10)	2150
Air inlet filters (pack/4)	2153A
Bacterial filter	1214

Please visit [www.philips.com/respironics](http://www.philips.com/respironics)



© 2011 Koninklijke Philips Electronics N.V.

All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

MCI 4103775 1080008 / bdc13041 \* APR 2011

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters  
 +33 1 47 28 30 82  
 Philips Respironics Asia Pacific  
 +65 6882 5282  
 Philips Respironics Australia  
 +61 (2) 9666 4444  
 Philips Respironics United Kingdom  
 +44 800 1300 845

Philips Respironics New Jersey, Inc.  
 5 Wood Hollow Road  
 Parsippany, NJ 07054 USA

Customer Service  
 +1 724 387 4000  
 800 345 6443 (toll free, US only)  
[www.philips.com/respironics](http://www.philips.com/respironics)